


# CORE SERIES

(P.N. 9539, Revision B3, November 2008)



**Easy Reference:**

Model name of the balance:	
Serial number of the unit:	
Software revision number (Displayed when power is first turned on):	
Date of Purchase:	
Name of the supplier and place:	

## 1.0 CONTENTS

<b>1.0</b>	<b>Contents</b> .....	<b>3</b>
<b>2.0</b>	<b>Introduction</b> .....	<b>3</b>
<b>3.0</b>	<b>Set Up</b> .....	<b>4</b>
3.1	Unpacking And Setting Up Your Balance .....	4
3.2	Installing Batteries .....	5
3.3	Overview - Display / Key Board.....	5
3.4	Locating And Protecting Your Balance .....	7
3.5	Storing When Not In Use .....	8
<b>4.0</b>	<b>Basic Operation</b> .....	<b>9</b>
4.1	Turning On The Balance .....	9
4.2	Zeroing / Tare .....	9
4.3	Weighing.....	10
4.4	Weighing Units .....	10
<b>5.0</b>	<b>Parameters</b> .....	<b>11</b>
5.1	Auto Power Off .....	11
5.2	Setting The Backlight .....	12
5.3	Enabling The Weighing Units .....	12
<b>6.0</b>	<b>Calibration</b> .....	<b>14</b>
<b>7.0</b>	<b>Trouble Shooting</b> .....	<b>15</b>
7.1	Error Messages.....	15
7.2	Replacement Parts And Accessories.....	16
<b>8.0</b>	<b>Specifications</b> .....	<b>17</b>
<b>9.0</b>	<b>Calibration Certificate</b> .....	<b>18</b>
<b>10.0</b>	<b>Warranty Statement</b> .....	<b>18</b>


## 2.0 INTRODUCTION

Thank you for purchasing your new Core™ Balance. Every Core Series balance has a range of features and options making them ideal for both the laboratory and classroom setting. We hope you enjoy using your new Core Balance.


## 3.0 SET UP

### 3.1 UNPACKING AND SETTING UP YOUR BALANCE

- 1) Remove the shipping protection screw and round cover as shown below using the Allen key supplied in the box. (It is advisable to keep the shipping protection screw and the round cover along with the packing box for using during shipping to another location in future).



- 2) Place the pan (which comes separately packed) in the receptacles on the top cover.



Use the shipping protection screw that you removed to screw the plastic pan into place. Be careful when putting the screw in that you do not use too much force as it can damaged the loadcell. You should just be able to use light force to secure the screw.

- 3) Place the stainless steel pan on top of the plastic pan.
- 4) Place the balance on a firm and flat surface for accurate weighing.

### 3.2 INSTALLING BATTERIES

Remove the battery compartment cover and insert six batteries (R6P/LR6 /AA size) into the battery compartment.

### 3.3 OVERVIEW - DISPLAY / KEY BOARD

