

## Core™ Compact Portable Balances



Speed, Performance, Value - The Right Balance



The Core series from Adam Equipment is the best choice for simple operation and economy. The tough durable ABS housing is designed to withstand field and laboratory environments, while being easy to clean and protected from accidental spills. With built-in ShockProtect™ protection Core balances can handle excessive overloads without a problem. Complete with a removable draft shield, AC adapter, plus integral security slot, Core balances are ready to work right out of the box.



### Technical Specifications


| Model                      | CQT 202          | CQT 251 | CQT 601 | CQT 1501 | CQT 2601         | CQT 2000 | CQT 5000 |
|----------------------------|------------------|---------|---------|----------|------------------|----------|----------|
| Maximum capacity           | 200g             | 250g    | 600g    | 1500g    | 2600g            | 2000g    | 5000g    |
| Readability                | 0.01g            | 0.1g    | 0.1g    | 0.1g     | 0.1g             | 1g       | 1g       |
| Repeatability (S.D.)       | 0.02g            | 0.1g    | 0.1g    | 0.1g     | 0.1g             | 1g       | 1g       |
| Linearity (±)              | 0.02g            | 0.1g    | 0.1g    | 0.2g     | 0.2g             | 1g       | 1g       |
| Typical stabilisation time | 3 Sec            | 2 Sec   |         | 3 Sec    |                  |          | 2 Sec    |
| Draught shield dimensions  | -                |         |         |          | -                |          |          |
| Pan size                   | 4.7" Ø / 120mm Ø |         |         |          | 5.7" Ø / 145mm Ø |          |          |

### General Specifications

|                            |  |
|----------------------------|--|
| Units of measurement       | Kilogram (kg), Gram (g), Pound (lb), Ounce (Oz), Carat (ct), Dram (d), Grain (GN), Newton (N), Taels T. (TL.T) |
| Power supply               | 12VAC 15mA   |
| Overall dimensions (wxdxh) | 6.8×10×3.4" / 173×255×86mm   |

### Key Accessories

- Security Lock and Cable
- In-use Cover
- Hard Carry Case with Lock



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)