



# Dune™ Compact Balances



Speed, Performance, Value - The Right Balance





### Technical Specifications

Model	DCT 201	DCT 2000	DCT 5000
Maximum capacity	200g	2000g	5000g
Readability	0.1g	1g	2g
Repeatability (S.D.)	0.1g	1g	2g
Linearity ( $\pm$ )	0.2g	2g	4g
Units of measurement	Gram (g), Gram (g), Ounce (Oz), Ounce (Oz), Pound (lb)	Gram (g), Ounce (Oz), Pound (lb)	
Pan size	5.7×5.7" / 145×145mm		

### General Specifications

Typical stabilisation time	2 Sec
Operating temperature	32°F - 104°F
Power supply	12V AC Adapter
Battery	6 x AA
Display	Large Display
Overall dimensions (wxdxh)	5.8×8.2×1.6" / 147×208×40mm



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)