


Highland™ Portable Precision Balances



Speed, Performance, Value - The Right Balance





Technical Specifications

Model	HCB 123	HCB 153	HCB 302	HCB 602H	HCB 602	HCB 1002	HCB 1502	HCB 3001
Maximum capacity	120g	150g	300g	600g	600g	1000g	1500g	3000g
Readability	0.001g	0.005g	0.01g	0.01g	0.02g	0.01g	0.05g	0.1g
Repeatability (S.D.)	0.003g	0.005g	0.01g	0.1g	0.02g	0.02g	0.05g	0.1g
Linearity (±)	0.003g	0.01g	0.02g	0.02g	0.02g	0.02g	0.1g	0.2g
Pan size	4.7" Ø / 120mm Ø							

General Specifications

Typical stabilisation time	2 Sec
Power supply	12VDC 800mA
Draught shield	Removable draught shield standard on all models
Draught shield dimensions	132mm
Overall dimensions (wxdxh)	6.8x9.9x3" / 172x251x75mm

Key Accessories

- Thermal Paper for Adam Printer (5 Pack)
- USB Cable
- Hard Carry Case with Lock
- Adam Printer
- In-use Cover
- RS-232 Cable
- Security Lock and Cable
- AdamDU - Data Collection Program



Adam Equipment - offices worldwide

www.adamequipment.com