

# Highland™ Portable Precision Balances



## Technical Specifications

| Model                | HCB 123          | HCB 153 | HCB 302 | HCB 602H | HCB 602 | HCB 1002 | HCB 1502 | HCB 3001 |
|----------------------|------------------|---------|---------|----------|---------|----------|----------|----------|
| Maximum capacity     | 120g             | 150g    | 300g    | 600g     | 600g    | 1000g    | 1500g    | 3000g    |
| Readability          | 0.001g           | 0.005g  | 0.01g   | 0.01g    | 0.02g   | 0.01g    | 0.05g    | 0.1g     |
| Repeatability (S.D.) | 0.003g           | 0.005g  | 0.01g   | 0.1g     | 0.02g   | 0.02g    | 0.05g    | 0.1g     |
| Linearity (±)        | 0.003g           | 0.01g   | 0.02g   | 0.02g    | 0.02g   | 0.02g    | 0.1g     | 0.2g     |
| Pan size             | 4.7" Ø / 120mm Ø |         |         |          |         |          |          |          |

## General Specifications

|                            |   |
|----------------------------|---|
| Typical stabilisation time | 2 Sec   |
| Power supply               | 12VDC 800mA                                     |
| Draught shield             | Removable draught shield standard on all models |
| Draught shield dimensions  | 132mm   |
| Overall dimensions (wxdxh) | 6.8×9.9×3" / 172×251×75mm                       |

## Key Accessories

- Thermal Paper for Adam Printer (5 Pack)
- USB Cable
- Hard Carry Case with Lock
- Adam Printer
- In-use Cover
- RS-232 Cable
- Security Lock and Cable
- AdamDU - Data Collection Program



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)