



North

Home

About Us

News

Applications

Events

Support

Career

English Choose your languageSearch

Applications

Moticam 480

Clinical > Digital > Matching Solutions



Moticam 480

Moticam 480

Multiple Output Microscopy Camera

Description

The Moticam 480 represents the ultimate in flexibility when it comes to choosing a microscopy camera solution.

Attach this camera to your microscope's eyepiece or to a trinocular view to live images on a PC, Projector or a TV Monitor simultaneously. The 480's triple output:

- USB Digital (640x480 pixels)

- S-Video (480 TV Lines)

- RCA (480 TV Lines)

Hardware-wired White Balance ensures accurate color representation for professional applications.

When connected to a computer, the Moticam 480 becomes an integrated Microscopy workstation with Motic's Images Plus 2.0 ML application.

Specifications

Camera & Accessories	
Chip Configuration	1/3" Multiple Output Digital / Analog Imaging Camera in portable casing
Data Output	Hi-Speed USB2.0 (640x480 pixels)
	S-Video (480 TV Lines)
	RCA (480 TV Lines)
Power Supply	Attachable Universal Power Supply 110V~240V
Included Software	Motic Images Plus 2.0ML for Windows Motic Images Plus 2.0 for OSX
Included Accessories	12mm Focusable Lens
	28, 30, 34 & 35mm Eyepiece Couplers
	Calibration Slide for accurate measurements

	Macro Tube S-Video, RCA & USB Cables
Min. Specifications	Windows: 2000/XP; 256MB RAM; P3; Built-in U Macintosh: OSX or higher; 256MB RAM; G3; Bu USB2.0* *Mac compatibility planned for middle of 2007. hardware update so existing units can not be u driver change.

Features

- Microscope Camera with 3 simultaneous outputs
- Hardware White Balancing for superios colour rendition
- Fits almost any microscope on the market either directly o
Eyepiece or onto the trinocular port (with relevant C-Moun
- All accessories and cables included

Asia Pacific | Europe | North America
Copyright © Motic®. All rights reserved.

www.motic.com