

## Technical Specifications


Model	PGL 203	PGL 2002	PGL 4001	PGL 6001	PGL 8001	PGL 10001	PGL 20001
Maximum capacity	200g	2000g	4000g	6000g	8000g	10kg	20kg
Readability	0.001g	0.01g	0.1g	0.1g	0.1g	0.1g	0.1g
Repeatability (S.D.)	0.002g	0.02g	0.1g	0.1g	0.1g	0.15g	0.1g
Linearity ( $\pm$ )	0.004g	0.04g	0.2g	0.2g	0.2g	0.4g	0.15g
Pan size	5.7×4.9" / 145×125mm	7.6×7.6" / 192×192mm				15.7×11.8" / 400×300mm	

## General Specifications

Typical stabilisation time	2-3 Sec
Power supply	AC adaptor and rechargeable battery pack
Calibration	External
Display	Backlit LCD display with dual digits and capacity tracker
Overall dimensions (wxdxh)	14.1×9.9×4.1" / 358×251×104mm

## Key Accessories

- USB to RS-232 Interface Cable
- Anti-Vibration Table
- In-use Cover
- Adam Printer
- AdamDU - Data Collection Program
- Security Lock and Cable
- Below Balance Hanger
- Dust Cover
- Thermal Paper for Adam Printer (5 Pack)
- RS-232 Cable



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)