

## ProBlot™ 35 Deluxe Rocking Platform

### Model 35 Large Rocking Platform



- Adjustable tilt and speed for exact level of agitation
- Safe for incubator and cold room use
- Durable construction
- Double platform option expands usable space

### Product Information

Labnet's ProBlot 35 Deluxe Rocker is made for those laboratories requiring a larger platform area and precise control. Both speed and degree of tilt are adjustable across a broad range to provide the exact level of agitation required for your application. It is equipped with a 2 hour timer that can also be set for continuous operation.

The heavy duty motor and sturdy rocking platform can support loads of up to 10 pounds. Closed loop, solid state speed control rocker operates at constant speed regardless of voltage or load fluctuations. Permanently lubricated ball bearings provide long life operation.

The ProBlot 35 Deluxe Rocker is available with a single platform, or for added workspace, two stacked platforms. It has a nonslip rubber mat to keep sample containers in place. The rocker can be used in environments from 4° to 65°C, making it ideal for incubator use.

### Back to Top

### Product Specifications

Speed range	4 - 160 rocks/minute
Timer	0 to 2 hours or continuous
Motion	Tilting
Tilt range	Adjustable, ±20°
Clearance between platforms	3.25 inches/8.25 cm
Maximum load	10 lb/4.6 kg
Ambient operating range	+4° to 65°C
Dimensions (W x D x H)	14.2 x 11.02 x 6.5 inches/36 x 28 x 16.5 cm
Weight	15 lb/6.7 kg
Electrical*	120V~, 60 Hz

\*Other voltages available upon request