

Centrifuges :: Minicentrifuges

Spectrafuge™ Mini Centrifuge

Join the mini revolution



- Supplied with both standard microtube and strip tube
- Ideal for quick spin downs and microfiltration
- Starts and stops in seconds
- Compact design

Product Information

This personal sized centrifuge is ideal for microfiltration and quick spin downs from the walls and caps of microcentrifuge tubes

Requiring less than six inches of bench space, the Spectrafuge Mini is supplied with rotors and adapters to accommodate 1.5 ml, 0.5 ml and 0.4 ml tubes as well as 0.2 ml strips and individual tubes. The standard and strip tube rotors are easily interchanged.

The easy-open, translucent lid utilizes a stainless steel hinge pin for durability. Five different lid colors are available including gray, blue, teal, purple and red.

An on/off switch is located on the side of the Spectrafuge Mini for starting and stopping operation. Alternatively, with the switch in the “on” position, the centrifuge can be started and stopped by closing and opening the lid.

Back to Top

Product Specifications

Maximum speed/RCF

6000rpm/2000 x g

Maximum capacity

Strip tube rotor	2 x 0.2 ml strips
Standard tube rotor	6 x 1.5/2.0 ml
Dimensions (W x D x H)	5 x 5 x 4.4 inches/12.7 x 12.7 x 11.2 cm
Weight	1 lb/0.45 kg
Electrical*	120V~, 60 Hz
Warranty	12 Months

1 Average max rcf in the strip tube rotor is 1177 x g *Other voltages available upon request