







[Products](#) [Contact Us](#) [Ordering](#) [OEM Offerings](#) [Labnet Main](#)
[Whats New](#) [Register Product](#) [Search](#) [Service](#) [Promotions](#)

[Request Quote](#) | [Specifications](#) | [Order Information](#)

Centrifuges :: Compact Laboratory Centrifuges

Z100A Clinical Centrifuge

[← Back to Category](#)

[← Back to Group](#)



- Max Capacity 6 x 15ml
- Variable Speed Up To 4,000 rpm (1,900 x g)
- Clear Lid With Safety Lock
- Adapters Availbale For All Common Blood Draw Tubes

Product Information

The Z100A features a 6 x 15 mL fixed angle rotor. The variable speed unit runs to 4,000 rpm/1,900 x g, suitable for a variety of applications in the clinical setting or research lab. The rotor accepts six 15mL or 10mL tubes directly. Optional adapters accept smaller tubes, including common blood draw tubes.

The Z100A is equipped with a clear lid for visualizing the rotor. A safety interlock cuts power to the centrifuge when the lid is opened. The rugged ABS housing adds durability to the centrifuge and is easy to clean.

[Back to Top](#)

Product Specifications

Speed Control	Variable Speed
Speed	1,000 - 4,000 rpm
G-force	118 - 1,900 x g
Timer	1-30min. or continuous
Dimensions	11 x 13.75 x 10.2 in
Dimensions	26.5 x 33 x 24.5 cm
Electrical	120V/50-60 Hz, 1 Amp

[Back to Top](#)

Ordering Information

Catalog Number	Item Description	Price
C0100-A	Clinical Z100-A Centrifuge With Rotor	Request Quote
C0230-17A	Adapters for 5ml (13 x 75mm) and 7ml (13 x 100mm)tubes, pack of 4	Request Quote
C0200-17A	Adapters for 5ml (13 x 75mm) and 7ml (13 x 100mm)tubes, pack of 6	Request Quote

[Back to Top](#)

Download related documents

[Z100A Instruction Manual](#)

North America | Toll Free: 888-LABNET1 (522-6381) | (732) 417-1750 | labnet@labnetlink.com

[Products](#) [What's New](#) [Contact Us](#) [Employment](#) [Register Product](#) [Ordering](#) [Search](#) [OEM Offerings](#) [Service](#) [Labnet Main](#) [Promotions](#)