

HERMLE
Labnet

COMPACT

CENTRIFUGE



- Advanced digital control with rpm/rcf conversion
- Ideal for all common tube sizes, 0.2 to 50 ml
- Useful for blood, urine and cell separations



HERMLE
Labnet



HERMLE
Z 206 A

HERMLE Z 206 A

The compact Z 206 A is useful for applications in both the research and clinical settings. Four rotor options allow this centrifuge to accept nearly every tube size from 0.2 ml to 50 ml. Speed is programmable to a maximum of 6,000 in rpm or 4,180 rcf, while acceleration and deceleration rates can be adjusted to 10 different levels. This rugged unit and accessories are built to last. For added durability, the 15 ml and 50 ml angle rotors are constructed of stainless steel.

Distinct Control Panel:

- Simple one handed operation
- Set and view speed in both rpm or g-force
- Pre-selection of rotor code for proper g-force conversion
- Control over all settings with one ergonomic adjustment knob
- Permanent indication of preset and actual values
- 10 acceleration and deceleration levels
- Pre-selection of running time from 10 sec to 99 h 59 min or cont.
- Splash proof digital display

Features:

- Microprocessor control with large LCD display
- Four available rotor options, capable of tube sizes up to 50ml
- Imbalance identification system with automatic shut-off
- Noise level under 60 dBA
- Stainless steel angle rotors for 15ml and 50ml conical tubes
- Compact, 11 x 14.6 x 10.2 in.
30 lb/15 kg

Rotors:



Item No: C0200-96
Capacity: 12 x 15 ml
Max. Rpm: 6,000 rpm
Max. Rcf: 4,180 x g
Radius: 10.4 cm
Angle: 28°C

Adapters:
C0200-17A
5ml 12 x 75mm
7ml 13 x 100mm
10ml 16 x 100mm



Item No: C0200-97
Capacity: 6 x 50 ml
Max. Rpm: 6,000 rpm
Max. RCF: 3,820 x g
Radius: 9.5 cm
Angle: 28°C

Adapters:
C0232-8B
15ml 17 x 125mm
10ml 16 x 100mm
C0232-9B
5ml 12 x 75mm
7ml 13 x 100mm



Item No: C0200-18
Capacity: 6 x 5 ml
Max. Rpm: 3,500 rpm
Max. RCF: 1,450 x g
Radius: 10.4 cm
Angle: 0 to 90°C



Item No: C0200-1
Capacity: 18 x 1.5 ml
Max. Rpm: 6,000 rpm
Max. RCF: 2,930 x g
Radius: 7.3 cm
Angle: 45°C

Adapters:
C1205
0.5ml
C1222
0.2ml
C1206
0.4ml

Ordering Information:

C0206-A Z206-A Compact Centrifuge

HERMLE
Labnet

Labnet
Labnet International, Inc.

PO Box 841 - Woodbridge, NJ 07095
<http://www.labnetlink.com>
Phone: 1-888-LABNET1 / 1-732-417-0700
Fax: 1-732-417-1750
Email: labnet@labnetlink.com

