



North

Home

About Us

News

Applications

Events

Support

Career

English Choose your languageSearch

Applications

Moticam 1000

Clinical > Digital > Matching Solutions



Moticam 1000

Moticam 1000

Affordable High-Resolution Live Imaging Camera

Description


The Moticam 1000 is a competitively priced high-resolution live in microscopy camera. With a live image of 1280x1024 pixels, this is projected in any lecture or classroom environment without fear of Enjoy crisp and clear live images.

Motic Images Plus is an imaging suite included for both PC and Mac platforms providing measurement, analysis and sharing tools for students and researchers alike.

With all necessary accessories included, the Moticam 1000 fits virtually any Microscope's eyepiece therefore allowing you to turn any microscope into a Digital Microscope.

Specifications

Camera & Accessories	
Chip Configuration	1/2" Live 1.3 Mega Pixel Imaging Chip in port
Data Output	USB2.0
Power Supply	Through USB Cable from Computer
Included Software	Motic Images Plus 2.0ML for Windows & Motic Images Plus 2.0 for Macintosh OSX
Included Accessories	16mm Focusable Lens 28, 30, 34 & 35mm Eyepiece Couplers Macro Tube Calibration Slide for accurate measurements Adaptor for Tripod
Min. Specifications	Windows: 2000 or XP; 256MB RAM; P4; Built-in Webcam Macintosh: OSX or higher; 256MB RAM; G4; iMac



Asia Pacific | Europe | North America
Copyright © Motic®. All rights reserved.

www.motic.com