


Moticam 2000/2300

- Both Moticam are compact and lightweight cameras that attach to virtually any microscope through the eyepiece, eyetube or onto a trinocular mount
- Easy plug & play USB2.0 Hi-Speed connection ensures that crisp and clear real-time images are displayed on the computer monitor turning any microscope into a Digital Microscope
- Use these cameras with your own software through a TWAIN interface or use the included Motic Images Plus 2.0 application software on both Windows and Macintosh systems
- The complete set of included accessories gives these hi-resolution cameras the ability to be used in many different applications from Macro to Micro viewing
- Live 2.0 and 3.0 Megapixel Resolution

Simple Plug & Play connection



Moticam 2000 &

Imaging Device

Effective Pixels

Max. Still Image