



North

Home

About Us

News

Applications

Events

Support

Career

English Choose your languageSearch

Applications

Clinical > Digital > Matching Solutions

Moticam 3000



3.3MP CCD

Professional grade Documentation Camera with optional Cooling I

Description

The Moticam 3000 series offers you a professional Documentation with a high-resolution CCD imaging device in a solid and heat-dis-housing.

As an option, the Moticam 3000C includes a Peltier Cooling device used for Fluorescence Microscopy Imaging.

The Moticam 3000 series includes our Motic Images Advanced 3.2 Suite that features RGB and Grey Scale Manual and Automatic Se as well as our Motic Images Assembly and Multi Focus Programs.

Moticam 3000

Specifications

Camera & Accessories	
Chip Configuration	1/2" Live 3.3 Megapixels Colour CCD Imaging portable casing
Data Output	IEEE 1394
Optional Cooling	Optional Peltier Cooling Mechanism to 10 Deg below ambient (Moticam 3000C)
Power Supply	Either through 6-pin to 6-pin IEEE Cable or 4-with external power supply
Included Software	Motic Images Advanced 3.2
Included Accessories	IEEE 1394 PCI Card Carrying/Storage Case
Min. Specifications	Windows: 2000 or XP; 256MB RAM; P4;