



TRUSTED PROTECTION®
INSIDE AND OUT

TECHNICAL DATA SHEET

ALOEFORM® X6 NITRILE EXAM GLOVES WITH ACTIVALOE®

PRODUCT DESCRIPTION

PRODUCT:	AloeForm® X6		
GRADE:	Exam		
MATERIAL:	Nitrile	Interior: Powder-Free	
THICKNESS:	Finger: 5.9mil / 0.15mm	Palm: 5.1mil / 0.13mm	Cuff: 4.0mil / 0.10mm
COLOR:	Light Green		
SURFACE:	Fully Textured		
LENGTH:	9.5" / 242mm	Weight: 5.5g	
FIT TYPE:	Ambidextrous		
CUFF TYPE:	Beaded		
WATERTIGHT AQL:	0.65		
ORIGIN COUNTRY:	China		


PHYSICAL PROPERTIES

	BEFORE AGING	AFTER AGING	ASTM METHOD
TENSILE STRENGTH:	20 MPa	20 MPa	D412-10
ELONGATION:	570%	550%	D412-10

ORDERING

	XS	S	M	L	XL	XXL
SIZE:	5.5 - 6.0	6.5 - 7.0	7.5 - 8.0	8.5 - 9.0	9.5 - 10	10.5 - 11.0
AVG. PALM WIDTH:	N/A / N/A	3.1"/80mm	3.4"/87mm	4.0"/101mm	4.5"/114mm	4.8"/123mm
PART NUMBER:	N129401	N129402	N129403	N129404	N129405	N129406
PACKAGING:	100 gloves per box, 10 boxes per case					

INTEGRATED TECHNOLOGY



- **MilTek™** offers durability and tactile sensitivity.
- **TouchTek™** protects skin and prevents contamination.
- Patented **ActivAloe™** soothes dry, chapped skin.

PERFORMANCE STANDARDS

EN374-1:2003



EN374-2



LEVEL 2

REGULATORY COMPLIANCE

EN 420: 2003
EN 455: 2000-2015
ASTM F1671-97
ASTM F739-12
REACH EU 1907/2006
DIRECTIVE 2011/65/EU
EU Reg 528/2012
CAT III COMPLEX DESIGN
FDA 21 CFR 177.2600
Component materials comply with
FDA regulations for food contact



EN455

CE 0493



EAC

TP TC 019/2011



Tested for use with Chemotherapy Drugs. Gloves used for protection against chemotherapy drug exposure should be selected specifically for the type of chemicals tested. For a complete list of the drugs and concentrations tested for, and testing method used, please call 510-429-8692.

RECOMMENDED APPLICATIONS

- Assembly
- Maintenance
- Fabrication
- Painting
- Cleaning
- Chemical
- Liquid Handling
- Oily or Wet Surfaces
- Laboratory
- Clinical

STORAGE

Keep out of sunlight. Store in a cool, dry place. Keep away from sources of ozone or ignition.

