



TRUSTED PROTECTION®
INSIDE AND OUT

TECHNICAL DATA SHEET

ALOESKIN®
NITRILE EXAM GLOVES WITH ACTIVALOE®

PRODUCT DESCRIPTION

PRODUCT:	AloeSkin®		
GRADE:	Exam		
MATERIAL:	Nitrile	Interior: Powder-Free	
THICKNESS:	Finger: 3.6mil / 0.09mm	Palm: 2.7mil / 0.07mm	Cuff: 2.3mil / 0.06mm
COLOR:	Light Green		
SURFACE:	Textured		
LENGTH:	9.5" / 242mm	Weight: 3.5g	
FIT TYPE:	Ambidextrous		
CUFF TYPE:	Beaded		
WATERTIGHT AQL:	1.5		
ORIGIN COUNTRY:	Malaysia		


PHYSICAL PROPERTIES

	BEFORE AGING	AFTER AGING	ASTM METHOD
TENSILE STRENGTH:	23 MPa	22 MPa	D412-10
ELONGATION:	520%	450%	D412-10

ORDERING

	XS	S	M	L	XL	XXL
SIZE:	5.5 - 6.0	6.5 - 7.0	7.5 - 8.0	8.5 - 9.0	9.5 - 10	10.5 - 11.0
AVG. PALM WIDTH:	N/A / N/A	3.4"/85mm	3.7"/96mm	4.1"/105mm	4.5"/114mm	N/A / N/A
PART NUMBER:	N193401	N193402	N193403	N193404	N193405	N/A
PACKAGING:	200 gloves per box, 10 boxes per case					

INTEGRATED TECHNOLOGY



- TouchTek™ protects skin and prevents contamination.
- Patented **ActivAloe™** soothes dry, chapped skin.

PERFORMANCE STANDARDS

EN374-1:2003



EN374-2



LEVEL 2

REGULATORY COMPLIANCE

EN 420: 2003
EN 455: 2000-2015
ASTM F1671-97
ASTM F739-12
REACH EU 1907/2006
DIRECTIVE 2011/65/EU
EU Reg 528/2012
CAT III COMPLEX DESIGN
FDA 21 CFR 177.2600
Component materials comply with
FDA regulations for food contact

CE
EN455

CE 0493

EAC
TP TC 019/2011

EAC

RECOMMENDED APPLICATIONS

- Assembly
- Maintenance
- Fabrication
- Painting
- Cleaning
- Chemical
- Liquid Handling
- Oily or Wet Surfaces
- Laboratory
- Clinical

STORAGE

Keep out of sunlight. Store in a cool, dry place.
Keep away from sources of ozone or ignition.

