
# **BioPette**™ PLUS

Labnet International, Inc.







- Dual volume adjustment with volume lock
- Universal tip compatibility
- Fully autoclavable



# **BioPette**™ PLUS

## Single Channel Pipette

Labnet's new BioPette Plus pipettes have been engineered to provide the highest level of quality, comfort, accuracy, and precision. Available in single, eight, and twelve-channel pipettes, they provide a range from 0.1µl to 10,000µl.

Each unit is supplied with a calibration certification and a three-year warranty. The BioPette Plus is completely autoclavable and UV resistant. The BioPette Plus chemical-resistant shafts have been designed for universal tip compatibility. The single channel has stainless steel tip ejectors which are easily removed and adjusted to fit any tip.





## Volume adjustment Finger rest using both pushbutton and thumbwheel Volume lock Ergonomic handle fits in both right Continuous and left hand volume adjustment Fully autoclavable 360° rotating and UV resistant manifold 20-200µI Individual spring loaded pistons Curved ejector bar Universal tip compatibility

# **BioPette**™ PLUS

### Multichannel Pipette

The BioPette Plus finger rest is ergonomically designed, allowing easy single-handed operation and fits comfortably in both right and left hands. The Biopette Plus also features a continuous volume adjustment, and a uniquely designed pushbutton and volume lock which prevents accidental volume change.

While having remarkably low plunger and volume change forces, Biopette Plus pipettes are highly accurate and precise. Each pipette is provided with calibration key for easy in-lab calibration.

BioPette Plus multichannel pipettes feature individually spring-loaded piston assemblies for positive tip sealing. They are provided with a curved ejector bar to produce low ejector forces, and a manifold that rotates 360° for ejecting and pipetting in any direction.







### Specifications and Ordering Information

| Catalog No. | Volume range                                                                                                                             | Channels | Accuracy      | Precision               | Tips      |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------|----------|---------------|-------------------------|-----------|
| SINGLE CHAN | INEL                                                                                                                                     |          |               |                         |           |
| PIPETTES*   |                                                                                                                                          |          |               |                         |           |
| 23940- 2    | 0.1-2µl                                                                                                                                  | 1        | ±12 to ±1.5%  | ≤6.0 to ≤0.7%           | 10µl      |
| P3940-10    | 0.5-10μl                                                                                                                                 | 1        | ±4.0 to ±0.5% | ≤4.0 to ≤0.4%           | 10μl      |
| P3940-20    | 2-20µl                                                                                                                                   | 1        | ±4.0 to ±0.8% | ≤3.0 to ≤0.4%           | 200µl     |
| P3940-100   | 10-100μl                                                                                                                                 | 1        | ±1.6 to ±0.8% | ≤0.8 to ≤0.3%           | 200µl     |
| P3940-200   | 20-200μΙ                                                                                                                                 | 1        | ±1.2 to ±0.6% | ≤0.6 to ≤0.3%           | 200µl     |
| P3940-1000  | 100-1000µl                                                                                                                               | 1        | ±1.0 to ±0.6% | ≤0.45 to ≤0.25%         | 1000µl    |
| P3940-5000  | 1000-5000µl                                                                                                                              | 1        | ±0.6 to ±0.5% | ≤0.3 to ≤0.2%           | 5000µl    |
| P3940-10000 | 1000-10000µl                                                                                                                             | 1        | ±2.5 to ±0.5% | ≤0.6 to ≤0.2%           | 10000µl   |
| MULTICHANN  | NEL                                                                                                                                      |          |               |                         |           |
| PIPETTES*   |                                                                                                                                          |          |               |                         |           |
| P4508-10    | 1-10µl                                                                                                                                   | 8        | ±8.0 to ±2.0% | ≤6.0 to ≤1.2%           | 10µl      |
| P4508-50    | 5-50µl                                                                                                                                   | 8        | ±4.0 to ±1.6% | ≤2.5 to ≤0.6%           | 200µl     |
| P4508-200   | 20-200μl                                                                                                                                 | 8        | ±3.0 to ±1.0% | ≤1.5 to ≤0.6%           | 200µl     |
| P4508-300   | 50-300µl                                                                                                                                 | 8        | ±1.6 to ±1.0% | ≤1.5 to ≤0.6%           | 200/300   |
| P4512-10    | 1-10µl                                                                                                                                   | 12       | ±8.0 to ±2.0% | ≤6.0 to ≤1.2%           | 10µl      |
| P4512-50    | 5-50µl                                                                                                                                   | 12       | ±4.0 to ±1.6% | ≤2.5 to ≤0.6%           | 200µl     |
| P4512-200   | 20-200μl                                                                                                                                 | 12       | ±3.0 to ±1.0% | ≤1.5 to ≤0.6%           | 200µl     |
| P4512-300   | 50-300µl                                                                                                                                 | 12       | ±1.6 to ±1.0% | ≤1.5 to ≤0.6%           | 200/300լ  |
| ACCESSORIES | 5                                                                                                                                        |          |               |                         |           |
| P3985       | Carousel stand for six single channel pipettes                                                                                           |          |               | BioPette Plus Four Pack |           |
| P3989       | Linear rack for six single or multichannel pipettes                                                                                      |          |               | BioPette Plus           | rour Pack |
| P5423       | Multichannel stand holds one 8 or 12 channel pipette                                                                                     |          |               |                         |           |
| P3940-SC    | Shelf clip holds one single or multichannel pipette                                                                                      |          |               |                         |           |
| P3986       | Shelf clip holds three single or multichannel pipettes                                                                                   |          |               |                         |           |
| P1000-10R   | Pipette Tips, 10μl rack of 96                                                                                                            |          |               |                         |           |
| P1000-10X-R | Pipette Tips, extended length, 10µl rack of 96                                                                                           |          |               |                         |           |
| P1000-200R  | Pipette Tips, 200µl rack of 96                                                                                                           |          |               |                         | 13        |
| P1000-300R  | Pipette Tips, 300µl rack of 96                                                                                                           |          |               |                         |           |
| P1000-1000R | Pipette Tips, 1000µl rack of 100                                                                                                         |          |               | 0                       |           |
| P1000-1200R | Pipette Tips, 1200µl rack of 96                                                                                                          |          |               |                         |           |
| P3998       | Pipette Tips, natural,                                                                                                                   |          | 200           |                         |           |
| K943        | Pipette Tips, clear, 10                                                                                                                  | A TOPE   |               |                         |           |
| P3940-SK4   | BioPette Plus Four Pack includes four pipettes:<br>0.5 to 10µl , 2.0 to 20µl , 20 to 200µl, 100 to 1000µl<br>Six position carousel stand |          |               | 188                     | t 4/      |
|             |                                                                                                                                          |          |               | 1111                    |           |
|             |                                                                                                                                          |          |               |                         |           |
|             | Three racks of Labnet BioFree™ Pipette Tips                                                                                              |          |               |                         |           |



P.O. Box 841 • Woodbridge, NJ 07095 • USA • website: http://www.labnetlink.com

U.S. Toll free: 888-LABNET1 • Fax: 732-417-1750 • Email: labnet@labnetlink.com

INTERNATIONAL Phone: +1-732-417-0700 • Fax: +1-732-417-1750

Email: international@labnetlink.com