



TRUSTED PROTECTION®
INSIDE AND OUT

TECHNICAL DATA SHEET

COLOR-Q® ORANGE NITRILE EXAM GLOVES WITH DUAL COLOR

PRODUCT DESCRIPTION

| | | | |
|-----------------|-------------------------|-----------------------|-----------------------|
| PRODUCT: | Color-Q® Orange | | |
| GRADE: | Exam | | |
| MATERIAL: | Nitrile | Interior: Powder-Free | |
| THICKNESS: | Finger: 5.1mil / 0.13mm | Palm: 4.3mil / 0.11mm | Cuff: 3.7mil / 0.09mm |
| COLOR: | Exterior: Orange | Interior: White | |
| SURFACE: | Fully Textured | | |
| LENGTH: | 9.5" / 242mm | Weight: 5.3g | |
| FIT TYPE: | Ambidextrous | | |
| CUFF TYPE: | Beaded | | |
| WATERTIGHT AQL: | 0.65 | | |
| ORIGIN COUNTRY: | China | | |


PHYSICAL PROPERTIES

| | BEFORE AGING | AFTER AGING | ASTM METHOD |
|-------------------|--------------|-------------|-------------|
| TENSILE STRENGTH: | 19 MPa | 18 MPa | D412-10 |
| ELONGATION: | 550% | 540% | D412-10 |

ORDERING

| | XS | S | M | L | XL | XXL |
|------------------|---------------------------------------|-----------|-----------|------------|------------|-------------|
| SIZE: | 5.5 - 6.0 | 6.5 - 7.0 | 7.5 - 8.0 | 8.5 - 9.0 | 9.5 - 10 | 10.5 - 11.0 |
| AVG. PALM WIDTH: | 2.9"/75mm | 3.4"/85mm | 3.7"/95mm | 4.2"/107mm | 4.5"/114mm | N/A / N/A |
| PART NUMBER: | N080671 | N080672 | N080673 | N080674 | N080675 | N/A |
| PACKAGING: | 100 gloves per box, 10 boxes per case | | | | | |

INTEGRATED TECHNOLOGY



- TouchTek™ protects skin and prevents contamination.
- **COLOR-Q®** enhances durability, engineered with two distinct polymer layers.

PERFORMANCE STANDARDS

EN374-1:2003



EN374-2



LEVEL 2

REGULATORY COMPLIANCE

EN 420: 2003
EN 455: 2000-2015
ASTM F1671-97
ASTM F739-12
REACH EU 1907/2006
DIRECTIVE 2011/65/EU
EU Reg 528/2012
CAT III COMPLEX DESIGN
UN NFPA 1999-2018
FDA 21 CFR 177.2600
Component materials comply with
FDA regulations for food contact

CE
EN455

CE 0493

EAC
TP TC 019/2011

EAC

RECOMMENDED APPLICATIONS

- Assembly
- Maintenance
- Fabrication
- Painting
- Cleaning
- Chemical
- Liquid Handling
- Oily or Wet Surfaces
- Laboratory
- Clinical

STORAGE

Keep out of sunlight. Store in a cool, dry place.
Keep away from sources of ozone or ignition.

