



DISTILLATION APPARATUS

DSA001




Figure 1

DESCRIPTION

The DSA001 Distillation Apparatus is a kit of the glassware and support materials needed for simple distillation experiments. It features high quality borosilicate glass components and robust stands and clamps. A storage box with a fitted foam insert is provided for the glassware. The components of the apparatus are shown in *Figure 1* and identified in the table below.

IDENTIFICATION OF COMPONENTS

- | | | | |
|----|--|-----|----------------------------------|
| 1. | Support Rod, 585 mm X 13 mm | 6. | Right Angle Clamp |
| 2. | 500ml Distillation Flask with Still Head | 7. | Extension Clamp with 1 Fixed Jaw |
| 3. | Recipient (250 ml Erlenmeyer Flask) | 8. | Support Base |
| 4. | Still Head Stopper | 9. | Extension Clamp, 2 Moveable Jaws |
| 5. | Graham-Type Condenser, 300 mm | 10. | Right Angle Clamp |