

Sample Tubes for Storage and Transportation



Globe Scientific's sample tubes have been designed for the storage and transport of biological material. They are produced from non-toxic polypropylene (PP) and offer a high degree of strength and clarity. The tubes feature external threads which

provide a smooth, uniform inner surface.

- Tubes withstand a temperature range of -196 °C to +121 °C
- Round bottom tubes can be centrifuged up to 17,000 x g
- Autoclavable (121 °C) in an upright position with caps loosened

Graduated with White Writing Area

Caps sold separately. See reverse for details.

(shown actual size)



Order Number	6050	6051	6052	6053	6054	6055	6058	6059
Volume	1.2 mL	2.0 mL	2.0 mL	3.0 mL	4.0 mL	4.0 mL	5.0 mL	5.0 mL
Height	40.9 mm	44.7 mm	47.4 mm	69.4 mm	72.5 mm	74 mm	88 mm	89 mm
Style	Self standing	Round bottom	Self standing	Self standing	Round bottom	Self standing	Round bottom	Self standing