



TRUSTED PROTECTION®
INSIDE AND OUT

TECHNICAL DATA SHEET

PH NATURAL®
pH BALANCING HAND HEALTH

PRODUCT DESCRIPTION

PRODUCT:	pH Natural®		
GRADE:	Exam		
MATERIAL:	Nitrile	POWDER: Powder-Free	
THICKNESS:	Finger: 4.5mil / 0.11mm	Palm: 3.6mil / 0.09mm	Cuff: 3.2mil / 0.08mm
COLOR:	Blue		
SURFACE:	Fully Textured		
LENGTH:	9.5" / 242mm	WEIGHT: 5.1g	
FIT TYPE:	Ambidextrous		
CUFF TYPE:	Beaded		
WATERTIGHT AQL:	0.65		
ORIGIN COUNTRY:	China		

PHYSICAL PROPERTIES

	BEFORE AGING	AFTER AGING	ASTM METHOD
TENSILE STRENGTH:	24 MPa	21 MPa	D412-10
ELONGATION:	580%	550%	D412-10

ORDERING

	XS	S	M	L	XL	XXL
SIZE:	5.5 - 6.0	6.5 - 7.0	7.5 - 8.0	8.5 - 9.0	9.5 - 10	10.5 - 11.0
AVG. PALM WIDTH:	3.0"/76mm	3.3"/85mm	3.7"/95mm	4.2"/107mm	4.5"/115mm	N/A/N/A
PART NUMBER:	N125331	N125332	N125333	N125334	N125335	N/A
PACKAGING:	100 gloves per box, 10 boxes per case					

INTEGRATED TECHNOLOGY



- **MilTek™** offers durability and tactile sensitivity.
- **TouchTek™** protects skin and prevents contamination.
- **pH Natural™** coating helps control the pH levels of the skin while wearing gloves.

PERFORMANCE STANDARDS

EN374-1:2003



EN374-2



LEVEL 2

REGULATORY COMPLIANCE

EN 420: 2003
EN 455: 2000-2015
ASTM F1671-97
ASTM F739-12
REACH EU 1907/2006
DIRECTIVE 2011/65/EU
EU Reg 528/2012
CAT III COMPLEX DESIGN
FDA 21 CFR 177.2600
Component materials comply with FDA regulations for food contact



TP TC 019/2011



RECOMMENDED APPLICATIONS

- Assembly
- Maintenance
- Fabrication
- Painting
- Cleaning
- Chemical
- Liquid Handling
- Oily or Wet Surfaces
- Laboratory
- Clinical

STORAGE

Keep out of sunlight. Store in a cool, dry place. Keep away from sources of ozone or ignition.

