

# Supercede® X8 Nitrile Exam Gloves



Supercede X8 powder-free nitrile exam gloves provide excellent chemical and biological barrier protection with an ultra-clean TouchTek™ surface.

Proprietary MilTek™ polymer technology delivers more uniform thickness providing unmatched tactile sensitivity at the fingers, high durability at the palm and cuff.

TouchTek™ surfaces reduce irritants and contaminants to protect your hands and sensitive items you may be handling.

Non-latex nitrile offers superior abrasion and puncture resistance and the fully textured grip allows better handling of small parts and slippery objects.

## Features and Benefits

- TouchTek™ Technology prevents skin irritation and contamination
- MilTek™ offers durability and tactile sensitivity
- Fully Textured Surface delivers excellent grip
- Pure Nitrile prevents latex allergic reactions
- Beaded Cuff allows easy donning and removal

## Typical Applications:

- Light Metal Fabricating
- Machinery Repair
- Cleaning & Maintenance
- Janitorial/Sanitation
- Painting

## Specifications:

Product Category:	Exam Gloves	
Material:	Nitrile	
Color:	Light Blue	
Inner Surface:	Smooth	
Outer Surface:	Textured	
Cuff Finish:	Beaded	
Fit Type:	Ambidextrous	
Length*:	240mm	
Finger Thickness:	0.18mm	7.09mil
Palm Thickness:	0.16mm	6.30mil
Cuff Thickness:	0.13mm	5.12mil
Weight*:	8.10g	
Online Processing:	Chlorination	
Tensile Strength Before Aging:	24 MPa	
Tensile Strength After Aging:	23 MPa	
Elongation Before Aging:	600%	
Elongation After Aging:	550%	
Watertight Barrier AQL:	1.5	
Powder Mass:	<2.0 mg/glove	
Sterile or Non-Sterile:	Non-Sterile	
Packaging Type:	50 Gloves/Dispenser Box	

\*Specification based on size medium gloves.

## Ordering Information:

Item No.	Size	Packaging Quantities
N098312	Small	50 Gloves/Box, 500/Case
N098313	Medium	50 Gloves/Box, 500/Case
N098314	Large	50 Gloves/Box, 500/Case
N098315	Extra Large	50 Gloves/Box, 500/Case
N098316	XX Large	50 Gloves/Box, 500/Case

## In Compliance With:

- EN 420 2003
- EN 455 2000-2009
- REACH EU 1907/2006
- Directive 2011/65/EU
- EU Reg 528/2012
- FDA 21 CFR 177.2600
- PPE CAT I 89/686/EEC  
For minimal risks only

## Meets Standards:

- ASTM D6319-10
- ASTM F1671-97
- ASTM F739-12

## Tested For Use With Chemotherapy Drugs:

Gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of drugs being used. Chemical permeation data of chemotherapy drugs and a Material Safety Data Sheet are available upon request. For chemical resistance data on chemotherapy drugs, please call +1.510.429.8692.



PROTECTIVE SINGLE-USE EXAMINATION GLOVES FOR EMERGENCY MEDICAL OPERATIONS IN ACCORDANCE WITH NFPA 1999, 2013 EDITION 3UM9. THIS GLOVE IS FOR SINGLE USE ONLY.

THIS GLOVE MEETS THE SINGLE-USE EMERGENCY MEDICAL EXAMINATION GLOVE REQUIREMENTS OF NFPA 1999, STANDARD ON PROTECTIVE CLOTHING AND ENSEMBLES FOR EMERGENCY MEDICAL OPERATIONS, 2013 EDITION. DO NOT REMOVE THIS LABEL.



NOT MADE  
WITH NATURAL  
RUBBER LATEX



SW Safety Solutions Inc. 33278 Central Ave. Unit 102 Union City, CA 94587 U.S.A. T: +1-510-429-8692 T: +1-510-487-5347

©2015 SW Safety Solutions Inc. SW and SUPERCEDÉ are registered trademarks of SW Safety Solutions Inc.

swsafety.com

M\_N09831X0815